

# **Paired Sequencing Reads**

**(Paired End and Mate Pair Sequencing)**


## Paired Sequencing Reads

**Paired Sequencing Reads** are created by sequencing DNA fragments from both ends.



**Paired Sequencing Reads** have been common from the early days of Sanger Sequencing.


# Paired Sequencing Reads



**Paired Sequencing Reads** must assemble predictable relative to each other.

# Paired Sequencing Reads


For “De Novo Assembly”, this enables combining Contigs to form Scaffolds.



**Paired Reads** that assemble in separate **Contigs** indicate **Order** and **Orientation**.

## Paired Sequencing Reads

With a **Reference Sequence** (“Resequencing” or “Mapping”),  
**Paired Reads** have a number of uses, including:




**Reference Sequence**

Assessing the accuracy of an assembly from the predictable assembly of the **Paired Reads**.

## Paired Sequencing Reads

With a **Reference Sequence** (“Resequencing” or “Mapping”),  
**Paired Reads** have a number of uses, including:




**Reference Sequence**

Detection of **Insertions** (assuming the **Template** size is fairly accurately known).

## Paired Sequencing Reads

With a **Reference Sequence** (“Resequencing” or “Mapping”),  
**Paired Reads** have a number of uses, including:



**Reference Sequence**

Detection of **Deletions** (assuming the **Template** size is fairly accurately known).

## Paired Sequencing Reads

With a **Reference Sequence** (“**Resequencing**” or “**Mapping**”),  
**Paired Reads** have a number of uses, including:



Resolving ambiguous assembly issues caused by repeat regions.

A single read might match several similar repeat regions.

If paired with a read outside the repeat, the correct match might be ascertained.



## Paired Sequencing Reads

Paired Sequencing Reads fall into two main categories.



### Paired End Sequencing:

Short template (~800 bp)

Reads point inwards from the ends of the template



### Mate Pair Sequencing:

Longer template (2,000 – 2,500 bp)

Reads point outwards from near the ends of the template

For further details, try this [excellent video](#)